

TECHNICAL DATA SHEET

Part Number	1001085
Description	Whatman Grade 1 Qualitative Filter Papers



Benefits :

- ✓ Whatman Grade 1 Qualitative Filter Papers are the most widely used filter paper for routine laboratory applications, qualitative analytical separations, and clarifying liquids. They are the standard grade cellulose filter papers for medium flow rates.
- ✓ Application : Used filter paper for routine applications with medium retention and flow rate. Covers a wide range of laboratory applications and is frequently used for clarifying liquids.
- ✓ Diameter : 85 mm
- ✓ Grade 1 Standard
- ✓ Material : Cellulose
- ✓ Nominal Ash Content : 0.06%
- ✓ Nominal air flow rate : 13 s/100 ml/in²
- ✓ Nominal basis weight : 87 g/m²
- ✓ Nominal thickness : 180 µm
- ✓ Typical Particle Retention in Liquid² : 11 µm
- ✓ Typical water flow rate : 57 ml/min
- ✓ Format : Circle

Sales unit :

100 Filters papers
