


TECHNICAL DATASHEET

Part number	<i>Consult sheet herabove</i>
Description	StarFrost slides With Poly-L-Lysine coating



Description :

- ✓ Made from soda lime glass
- ✓ Cleaned, unfated, polished
- ✓ Poly-L-Lysine coating on one side of the slide
- ✓ Improved and secured tissular adhesion
- ✓ Meet standard ISO 8037/1

- ✓ Dimensions : 26 x 76 mm
- ✓ Thickness : 1,0 mm
- ✓ Marking area : 20 mm

Color	White	Yellow	Blue	Pink	Green	Orange
Coupé	100319	100320	100323	100321	100322	100324
Ground 90°	100331	100332	100335	100333	100334	100366
Ground 45°	100325	100326	100329	100327	100328	100330
Bevelled	100343	100344	100347	100345	100346	100348

Sale unit :
50 slides