

TECHNICAL DATA SHEET

Part Number	104433
Description	Microcentrifuge tube 2 ml



Benefits :

- ✓ Made from polypropylene
- ✓ Screw cap with syringe piecable-port
- ✓ Skirted bottom
- ✓ Molded-in graduation
- ✓ Colour : natural

Certificate :

- ✓ Non sterile
- ✓ Autoclavable

Technical specifications :

- ✓ Volume : 2 ml
- ✓ Centrifugable up to 17 000 g
- ✓ Temperature range : -196°C to +121°C

Sales unit :

500 Microcentrifuge tubes 2 ml
