

Migration guide

New names

Nitrocellulose membranes

Western Protran BA85	→	Amersham Protran NC 0.45
Western Protran BA88	→	Amersham Protran NC 0.2
Western Protran BA79	→	Amersham Protran NC 0.1
Amersham Hybond ECL 0.45 µm	→	Amersham Protran Premium NC 0.45
Amersham Hybond ECL 0.2 µm	→	Amersham Protran Premium NC 0.2
Western OptiMem SA035	→	Amersham Protran Supported NC 0.45
Western OptiMem BA83	→	Amersham Protran Supported NC 0.2

PVDF membranes

Amersham Hybond P	→	Amersham Hybond PVDF 0.45
Amersham Hybond P40	→	Amersham Hybond PVDF 0.2
Amersham Hybond LFP	→	Amersham Hybond LFP PVDF 0.2
Western Western S	→	Amersham Hybond SE0 PVDF 0.2

With Amersham Western blotting membranes, GE Healthcare Life Sciences continue to provide scientists with excellent membranes with an expanded selection of formats, rolls, precuts, and sandwiches.

For more information visit:
www.gewesternblotting.com/western-blotting-product

Also find out about our extensive range of imagers:
www.gewesternblotting.com/western-blotting-product/imagers

Or contact your local representative

GE Healthcare Life Sciences is a member of the GE Healthcare Company.
Amersham, Hybond, Protran, and Western are trademarks of GE Healthcare Company.
© 2012 GE Healthcare Company. All rights reserved.
GE Healthcare Life Sciences is a registered trademark of GE Healthcare Company.
GE Healthcare Life Sciences is a registered trademark of GE Healthcare Company.

GE Healthcare Life Sciences AB, Nydalen 30, 751 84 Uppsala, Sweden

GE Healthcare
Life Sciences

Two great brands
join forces

New Amersham™ Western
blotting membranes



imagination at work

2M-0076-01-04 10/2012

New Amersham Western blotting membranes

by GE Healthcare Life Sciences


Two great brands
join forces



- Optimized for chemiluminescent and fluorescent detection
- Excellent protein binding capacity over a wide size range
- Get 30% more with the new pack sizes
- Physical and chemical composition remains the same

The trusted Hybond and Whatman brands combine to give Amersham Western blotting membranes, offering sustained high quality and performance.

New



Amersham membrane sandwiches—designed for convenience and ease-of-use

New



Amersham Hybond PVDF 0.2—specifically optimized for difficult Western blotting applications