

Technical Datasheet

Art.-No.: 1120

Buffer solution pH 10.00 (pH 10.00 (±0.02/20 °C))
traceable to SRM from NIST

For laboratory use.

Parameter	Value
Density	1.01 g/ml
H ₃ BO ₃ -KCl-NaOH	
pH (0 °C)	10.24
pH (10 °C)	10.10
pH (20 °C)	10.00
pH (25 °C)	9.95
pH (30 °C)	9.90
pH (40 °C)	9.82
pH (50 °C)	9.75
pH (60 °C)	9.68
pH (70 °C)	9.62
pH (80 °C)	9.55
pH (90 °C)	9.49