

**1. Part Numbers:** 8131A-50

**2. Photograph:**



**3. Product Information:**

Description	Sample Containers 50ml
Materials of Construction:	Body Polypropylene (PP) Lid High Density Polyethylene (HDPE)
Manufacturing Environment:	Controlled environment within factory
Nominal Volume	50ml
Brimful Volume	70ml
Colour:	Translucent with red lid
Label Area:	50 x 25mm
Graduation Marks:	None
Overall Height (with lid on):	59mm
Internal Neck Diameter:	43mm
Internal Base Diameter:	39mm
Body Height:	55mm
Nominal Weight:	12g
Individually Bagged?	No
Pre-sterilised?	No
Number of Containers per Box:	384
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 1-Mar-22