


	<b>Product Specification</b>	Specification Number	8131A-100
		Revision	01

Part No. 8131A-100  
Description Sample Containers 100ml



Nominal Volume 100ml  
Brimful Volume 145ml  
Colour: Translucent with red lid  
Label Area: N/A  
Graduation Marks: Every 20ml  
Material of Construction: Body Polypropylene (PP)  
Lid High Density Polyethylene (HDPE)  
Overall Height (with lid on): 76mm  
Internal Neck Diameter: 56mm  
Body Height: 72mm  
Nominal Weight: 14g  
Individually Bagged? No  
Pre-sterilised? No  
Number of Containers per Box: 280  
Recommended Storage Conditions: Dry and ambient temperature

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 14-Feb-14
----------------------------	--	--