

**1. Part Numbers:** 8126H-1000A**2. Photograph:****3. Product Information:**

Description	Wide Mouth HDPE Bottle 1000ml Amber
Materials of Construction:	Bottle HDPE (conforms to FDA 177.1520 & USP Class VI) Cap PP (conforms to FDA 177.1520 & USP Class VI)
Manufacturing Environment:	Controlled environment within factory
Nominal Volume:	1000ml
Brimful Volume:	1075ml
Nominal Height:	201mm
Diameter:	92mm
Internal Neck Diameter:	52mm
Weight:	102g
Pack Size:	6 bottles per pack
Method of Sterilisation:	N/A (these bottles are non-sterile)
BSE/TSE Free:	Yes
Supplied with lid screwed on:	Yes
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:

Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Oct-20

**Sampling Systems Ltd.**

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.

Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com