

1. Part Numbers: 8132P-10T

2. Photograph:



3. Product Information:

Description	Carboy 10L with Stopcock		
Material of Construction:	10L Carboy	PP	(conforms to FDA 177.1520 & USP Class VI)
	Cap	PP	(conforms to FDA 177.1520 & USP Class VI)
	Ring Seal	TPE	(conforms to FDA 177.1810 & USP Class VI)
	Stopcock	PP	(conforms to FDA 177.1520 & USP Class VI)
Nominal Volume:	10 Litres		
Nominal Height:	390mm		
Max. Diameter:	250mm		
Internal Neck Diameter:	65mm		
Weight:	1.3 kg		
Graduated:	Yes		
Autoclavable:	Yes		
Pack Size:	1		
Recommended Storage Conditions:	Dry and ambient temperature		

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 25-Jan-21