

**1. Part Numbers:** 8132P-20T

**2. Photograph:**



**3. Product Information:**

Description	Carboy 20L with Stopcock		
Material of Construction:	20L Carboy	PP	(conforms to FDA 177.1520 & USP Class VI)
	Cap	PP	(conforms to FDA 177.1520 & USP Class VI)
	Ring Seal	TPE	(conforms to FDA 177.1810 & USP Class VI)
	Stopcock	PP	(conforms to FDA 177.1520 & USP Class VI)
Nominal Volume:	20 Litres		
Nominal Height:	530mm		
Max. Diameter:	290mm		
Internal Neck Diameter:	65mm		
Weight:	2.0 kg		
Graduated:	Yes		
Autoclavable:	Yes		
Pack Size:	1		
Recommended Storage Conditions:	Dry and ambient temperature		

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 25-Jan-21



**Sampling Systems Ltd.**

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.

Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com