

	<h1>Product Specification</h1>	Specification Number	8111H-200S
		Revision	01

Part No. 8111H-200S

Description SteriWare Powder Funnel – 200mm diameter

Material of Construction: High Density Polyethylene (HDPE) - virgin

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004



Method of Construction: Single piece, injection moulded

Max. Diameter: 200mm

Overall Height: 165mm

OD at Bottom of Spout: 43mm

ID at Bottom of Spout: 40mm

Nominal Weight: 86g

Moulding & Packing Environment: Class 100,000 Cleanroom

Individually Bagged? Yes (heat sealed PE bag)

Method of Sterilisation: Gamma irradiated sterile (25 to 45kGy)


Material BSE/TSE Free Yes

Number of Funnels per Box: 10

BSE/TSE Free: Yes

Recommended Storage Conditions: Dry and ambient temperature

SteriWare is a registered trade marks belonging to Sampling Systems Ltd.

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 3-Mar-16
----------------------------	--	---