

**1. Part Number:** 8064H-500S

**2. Photograph:**



### 3. Product Information:

Description:	SteriWare® Powder Probe 500mm - Sterile
Use:	Single use powder sampler
Material of Construction:	High Density Polyethylene (HDPE) – virgin <ul style="list-style-type: none"><li>• Conforms to FDA 21 CFR 177.1520</li><li>• Conforms to EU Regulations 10/2011</li><li>• Conforms to EC Regulations 1935/2004</li><li>• BSE/TSE free</li></ul>
Method of Construction:	Ultrasonic welding
Nominal Length:	615mm
Blade Length:	497mm
Diameter of Blade:	22mm
Width of opening in blade:	10mm
Nominal Weight of Sampler:	89g
Moulding & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Quantity per Box:	100
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	Suitable for use in ATEX rated areas – Contact Sampling Systems for a copy of the report to confirm suitability for your application

### 4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable sampler. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 12-Mar-20

DOMINIQUE DUTSCHER SAS



**Sampling Systems Ltd.**

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.

Tel: +44(0)1675 466992 Email: [info@sampling.com](mailto:info@sampling.com) Web: [www.sampling.com](http://www.sampling.com)