

1. Part Numbers: 5035S-1000

2. Photograph:



3. Product Information:

Description	Liquid Master 1000mm – Stainless Steel	
Materials of Construction:	Tip	316L Stainless Steel
	Outer Tube	316L Stainless Steel
	Inner Rod	316L Stainless Steel
	Piston	316L Stainless Steel
	Handle	316L Stainless Steel
	Seal	Silicone (conforms to FDA CFR 177.2600)
	Spring	Stainless Steel
Method of Construction:	All components screw together	
External Surface Finish:	Better than 0.5 microns Ra	
Nominal Length:	1000mm	
Nominal Weight:	2.3kg	
Diameter of Outer Tube:	25mm	
Diameter of Hole in Tip:	6mm	
Recommended Storage Conditions:	Dry and ambient temperature	
Cleaning:	This product should be cleaned before first used	

4. Approval:

Signed:

Name: Kuldeep Flora

Position: Quality Manager


Date: 13-Oct-20



Sampling Systems Ltd.

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.

Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com



DOMINIQUE DUTSCHER SAS



ISO
9001 : 2015
REGISTERED
Cert No. 14128455

Sampling Systems Ltd.

4-5 Forge Mills Park, Coleshill. Warwickshire. B46 1JH. UK.

Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com