

1. Part Numbers: 5035S-1500

2. Photograph:



3. Product Information:

Description	Liquid Master 1500mm – Stainless Steel														
Materials of Construction:	<table border="0"> <tr><td>Tip</td><td>316L Stainless Steel</td></tr> <tr><td>Outer Tube</td><td>316L Stainless Steel</td></tr> <tr><td>Inner Rod</td><td>316L Stainless Steel</td></tr> <tr><td>Piston</td><td>316L Stainless Steel</td></tr> <tr><td>Handle</td><td>316L Stainless Steel</td></tr> <tr><td>Seal</td><td>Silicone (conforms to FDA CFR 177.2600)</td></tr> <tr><td>Spring</td><td>Stainless Steel</td></tr> </table>	Tip	316L Stainless Steel	Outer Tube	316L Stainless Steel	Inner Rod	316L Stainless Steel	Piston	316L Stainless Steel	Handle	316L Stainless Steel	Seal	Silicone (conforms to FDA CFR 177.2600)	Spring	Stainless Steel
Tip	316L Stainless Steel														
Outer Tube	316L Stainless Steel														
Inner Rod	316L Stainless Steel														
Piston	316L Stainless Steel														
Handle	316L Stainless Steel														
Seal	Silicone (conforms to FDA CFR 177.2600)														
Spring	Stainless Steel														
Method of Construction:	All components screw together														
External Surface Finish:	Better than 0.5 microns Ra														
Nominal Length:	1500mm														
Nominal Weight:	2.8kg														
Diameter of Outer Tube:	25mm														
Diameter of Hole in Tip:	6mm														
Recommended Storage Conditions:	Dry and ambient temperature														
Cleaning:	This product should be cleaned before first used														

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Oct-20