

1. Part Number: 8220A-003

2. Label Design:



3. Product Information:

| | |
|---------------------------------|---|
| Description: | Quality Control Label - Yellow - 150 x 150mm |
| Adhesive: | Ultra Adhesive - High strength foil backed paper labels |
| Colours: | Yellow printing on a white label |
| FDA Acceptable Adhesive? | No |
| Dimensions: | 150 x 150mm |
| Number of Labels per Roll: | 250 (approx). |
| Weight: | 1300g per roll |
| Method of Sterilisation: | N/A |
| Recommended Storage Conditions: | Dry and ambient temperature |

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 23-Sep-19



Sampling Systems Ltd.

4-5 Forge Mills Park, Coleshill, Warwickshire. B46 1JH. UK.
Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com