

**1. Part Number:** 8220A-32

**2. Label Design:**



The image shows a rectangular orange label with white text and input fields. At the top, the word "QUARANTINE" is printed in large white letters. Below it, there is a field for "Container" followed by a small square box and the word "of" followed by another small square box. The label is divided into several sections with white input boxes: "Material", "Batch / Lot Number", "Expiry Date", "Part Number", "Quantity", and "Initials & Date". On the left side, there is vertical text: "Sampling Systems Ltd. Tel +44 (0)1675 466992 E-mail: info@sampling.co.uk". On the right side, there is vertical text: "Re-order code 8220A-32 Visit our website at: www.sampling.co.uk".

### 3. Product Information:

**Description:** 'Quarantine' Label - 95 x 95mm  
**Adhesive:** Ultra Adhesive - High strength foil backed paper labels  
**Colours:** Orange printing on a white label  
**FDA Acceptable Adhesive?** No  
**Dimensions:** 95 x 95mm  
**Number of Labels per Roll:** 500 (approx).  
**Weight:** 1010g per roll  
**Method of Sterilisation:** N/A  
**Recommended Storage Conditions:** Dry and ambient temperature

### 4. Approval:

Signed:



**Name:** Kuldeep Flora

**Position:** Quality Manager

**Date:** 23-Sep-19