

	<h1 style="margin: 0;">Sampler Specification</h1>	Specification Number	8085A-100S
		Revision	05

Part No. 8085A-100S

Description Disposable Sack Sampler – 400mm – Sterile

Material of Construction: High Density Polyethylene (HDPE) - virgin

- Conforms to FDA CFR 177.1520
- Conforms to EU Regs 10/2011
- Conforms to EC Regs 1935/2004



Nominal Length: 400mm

Diameter of Outer Tube: 25mm

Nominal Weight of Sampler: 53g


Individually Bagged? Yes (heat sealed PE bag)

Method of Sterilisation: Gamma Irradiation (25 to 45 KGy)

Number of Samplers per Box: 100

BSE/TSE Free: Yes

Recommended Storage Conditions: Dry and ambient temperature

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 25-Nov-15
Page 1 of 1	Sampling Systems Ltd. 4-5 Forge Mills Park, Coleshill, Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Fax: +44(0)1675 466994 Email: info@sampling.com Web: www.sampling.com	