

**1. Part Number:** 8080A-900

**2. Product Information:**

Description: PureBag 900ml - Sterile

Use: Single use, disposable sampling bags

Material of Construction: Light Density Polyethylene (LDPE) - virgin

- Conforms to FDA CFR 177.1520
- BSE/TSE free

Closure Type: 2 steel wires (non product contact)

- Covered with FDA approved pressure sensitive tape with a backing of polypropylene or PVC

Bag Dimensions: 114 x 305mm

Thickness of Bag: 62 microns

Nominal Bag Volume: 900ml

Method of Sterilisation: Sterility is assured due to extrusion temperatures exceeding 220°C which eliminates any present bacteria. During manufacturing all bags are sealed at both top and bottom in a continuous format before being cut into individual bags. This method assures that the bag's interior will never come in contact with its surroundings. Periodic validation tests, based on USP 34 – Sterility testing, are performed in the manufacturing facility to guaranty the sterility of the processes.

Packaging: 1000 bags per case (2 x 500)

Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 5 years from the date of manufacture

ATEX: This product has not been ATEX tested

**3. Country of Origin:** Made and packed in Canada.

**4. Repeated Use:** This is a single-use, disposable product and should not be re-used.

**5. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 7-Oct-20

Note: The information given in this document has been given to us by our suppliers and is passed on in good faith, we believe the information to be correct.