

1. Part Number: 8069H-1000

2. Photograph:



3. Product Information:

Description: SteriWare® Mono Slot Sampler – 1000mm
Use: Single use, disposable powder sampler
Material of Construction: HDPE (High Density Polyethylene) - virgin

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

Method of Construction: Ultrasonic welding
Nominal Length: 968mm
Diameter of Blade: 22mm
Diameter of Outer Tube: 26mm
Number of Sampling Slots: 1
Slot Positions:



Sampling Slot Width: 8mm
Sampling Slot Length: 60mm
Nominal Weight of Sampler: 240g
Moulding Environment: Class 100,000 Cleanroom
Welding & Primary Packing Environment: Class 10,000 Cleanroom
Individually Bagged? Yes (heat sealed PE/PET bag)
Method of Sterilisation: N/A
Quantity per Box: 20
Recommended Storage Conditions: Dry and ambient temperature
Shelf Life: 5 years from the date of manufacture
ATEX: This product has not been ATEX tested



4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 22-Nov-22