

1. Part Number: 8073H-500

2. Photograph:



3. Product Information:

4. Description: SteriWare Slot Sampler - 500mm

5. Use: Single use, disposable powder sampler

Material of Construction: High Density Polyethylene (HDPE) - virgin

- Conforms to FDA 21 CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

Method of Construction: Ultrasonic welding


Overall Length: 620mm

Diameter of Outer Tube: 26mm

Diameter of Inner Blade: 22mm

Nominal Weight of Sampler: 148g

Slot Positions:



Number of Sampling Slots: 3

Sampling Slot Width: 8mm

Assembly & Primary Packing Environment: Class 100,000 Cleanroom

Individually Bagged? Yes (heat sealed PE/PET bag)

Method of Sterilisation: N/A

Number of Samplers per Box: 20

Recommended Storage Conditions:

Dry and ambient temperature

Shelf Life:

5 years from the date of manufacture

ATEX:

Suitable for use in ATEX rated areas – Contact Sampling Systems for a copy of the report to confirm suitability for your application

6. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

change.control@sampling.com**7. Country of Origin:**

Made and packed in the UK.

8. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

9. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Mar-20

10.