


**1. Part Number:** 8073H-500S

**2. Photograph:**



**3. Product Information:**

4. Description: SteriWare Slot Sampler - 500mm - Sterile  
 5. Use: Single use, disposable powder sampler  
 Material of Construction: High Density Polyethylene (HDPE) - virgin  
 • Conforms to FDA 21 CFR 177.1520  
 • Conforms to EU Regulations 10/2011  
 • Conforms to EC Regulations 1935/2004  
 • BSE/TSE free  
 Method of Construction: Ultrasonic welding  
 Overall Length: 620mm  
 Diameter of Outer Tube: 26mm  
 Diameter of Inner Blade: 22mm  
 Nominal Weight of Sampler: 148g  
 Slot Positions:



Number of Sampling Slots: 3  
 Sampling Slot Width: 8mm  
 Assembly & Primary Packing Environment: Class 100,000 Cleanroom  
 Individually Bagged? Yes (heat sealed PE/PET bag)  
 Method of Sterilisation: Gamma Irradiation (25-45kGy)  
 Number of Samplers per Box: 20

Recommended Storage Conditions:

Dry and ambient temperature

Shelf Life:

5 years from the date of manufacture

ATEX:

Suitable for use in ATEX rated areas – Contact Sampling Systems for a copy of the report to confirm suitability for your application

**6. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

[change.control@sampling.com](mailto:change.control@sampling.com)**7. Country of Origin:**

Made and packed in the UK.

**8. Repeated Use:**

This is a single-use, disposable sampler. Repeated use will lead to failure.

**9. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Mar-20