

**1. Part Number:** 8073H-1000S

**2. Photograph:**



**3. Product Information:**

4. Description: SteriWare Slot Sampler - 1000mm - Sterile

5. Use: Single use, disposable powder sampler

Material of Construction: High Density Polyethylene (HDPE) - virgin

- Conforms to FDA 21 CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

Method of Construction: Ultrasonic welding


Overall Length: 968mm

Diameter of Outer Tube: 26mm

Diameter of Inner Blade: 22mm

Nominal Weight of Sampler: 234g

Slot Positions:



Number of Sampling Slots: 3

Sampling Slot Width: 8mm

Assembly & Primary Packing Environment: Class 100,000 Cleanroom

Individually Bagged? Yes (heat sealed PE/PET bag)

Method of Sterilisation: Gamma Irradiation (25-50 kGy)

Number of Samplers per Box: 20

Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 5 years from the date of manufacture

ATEX: Suitable for use in ATEX rated areas – Contact Sampling Systems for a copy of the report to confirm suitability for your application



**Sampling Systems Ltd.**

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.

Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com

**6. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**7. Country of Origin:**

Made and packed in the UK.

**8. Repeated Use:**

This is a single-use, disposable sampler. Repeated use will lead to failure.

**9. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Mar-20