

1. Part Number: 8073J-1000

2. Photograph:



3. Product Information:

Description:	SteriWare Jumbo Slot Sampler - 1000mm
Use:	Single use, disposable powder sampler
Material of Construction:	High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> • Conforms to FDA 21 CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Ultrasonic welding
Overall Length:	986mm
Diameter of Outer Tube:	40mm
Diameter of Inner Tube:	35mm
Nominal Weight of Sampler:	297g
Number of Sampling Slots:	3
Sampling Slot Width:	8mm
Sampling Slot Length:	60mm
Assembly & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A
Number of Samplers per Box:	20
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	The device has not ben ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement. To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Mar-20

DOMINIQUE DUTSCHER SAS