

**1. Part Number:** 8050H-1500

**2. Photograph:**



### 3. Product Information:

**Description:**

SteriWare LiquiThief - 1500mm

**Use:**

Single use, disposable liquid sampler

**Material of Construction:**

High Density Polyethylene (HDPE) - virgin

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

O-ring Seal – virgin Silicone Rubber

- Conforms to FDA CFR 177.2600

**Method of Construction:**

Ultrasonic welding

**Nominal Length**

Fully closed:

1495mm

Fully extended:

2830mm

**Diameter of Outer Tube:**

21mm

**Diameter of hole in tip:**

1.5mm

**Width of Handle:**

103mm

**Nominal Weight of Sampler:**

240g

**Assembly & Primary Packing Environment:** Class 100,000 Cleanroom

**Individually Bagged?**

Yes (heat sealed PE/PET bag)

**Method of Sterilisation:**

N/A

**Quantity per Box:**

20

**Recommended Storage Conditions:**

Dry and ambient temperature

**Shelf Life:**

5 years from the date of manufacture

**ATEX:**

This device has not been ATEX tested

### 4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

[change.control@sampling.com](mailto:change.control@sampling.com)



**Sampling Systems Ltd.**

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.

Tel: +44(0)1675 466992 Email: [info@sampling.com](mailto:info@sampling.com) Web: [www.sampling.com](http://www.sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable sampler. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 6-Oct-20