


1. Part Number: 8055H-500

2. Photograph:



3. Product Information:

Description:

SteriWare® ViscoThief - 500mm

Use:

Single use, disposable liquid sampler (for high viscosity liquids)

Material of Construction:

High Density Polyethylene (HDPE) - virgin

- Conforms to FDA 21 CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- Conforms to REACH regulation 1907/2006
- BSE/TSE free

O-ring Seal – virgin Silicone Rubber

- Conforms to FDA CFR 177.2600

Method of Construction:

Ultrasonic welding

Nominal Length

Fully closed:

523mm

Fully extended:

881mm

Diameter of Outer Tube:

21mm

Diameter of hole in tip:

10mm

Width of Handle:

103mm

Nominal Weight of Sampler:

83g

Assembly & Primary Packing Environment: Class 100,000 Cleanroom

Individually Bagged?

Yes (heat sealed PE/PET bag)

Method of Sterilisation:

N/A

Number of Samplers per Box:

20

Recommended Storage Conditions:

Dry and ambient temperature

Shelf Life:

5 years from the date of manufacture

ATEX:

Suitable for use in ATEX rated areas – Contact Sampling Systems for a copy of the report to confirm suitability for your application



Sampling Systems Ltd.

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.

Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 6-Oct-20