

**1. Part Number:** 8051H-450S

**2. Photograph:**



**3. Product Information:**

**Description:**

SteriWare® Jumbo Pipette 450mm - sterile

**Use:**

Single use, disposable liquid sampler

**Material of Construction:**

High Density Polyethylene (HDPE) - virgin

- Conforms to FDA 21 CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

A blue dye is added to one of the tips. The blue dye is certified:

- FDA acceptable for food contact
- Free of materials of animal origin

**Method of Construction:**

Ultrasonic welding

**Nominal Length**

420mm

**Diameter of Outer Tube:**

21mm

**Diameter of hole in tip:**

Handle End (Blue) = 1.5mm  
Product End (Natural) = 3mm

**Nominal Weight of Sampler:**

50g

**Assembly & Primary Packing Environment:** Class 100,000 Cleanroom

**Individually Bagged?** Yes (heat sealed PE/PET bag)

**Method of Sterilisation:** Gamma Irradiation (25 kGy)

**Number of Samplers per Box:** 50

**Recommended Storage Conditions:** Dry and ambient temperature

**Shelf Life:** 5 years from the date of manufacture

**ATEX:** Suitable for use in ATEX rated areas – Contact Sampling Systems for a copy of the report to confirm suitability for your application

**Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**4. Country of Origin:**

Made and packed in the UK.

**5. Repeated Use:**

This is a single-use, disposable sampler. Repeated use will lead to failure.

**6. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 10-May-19