

**1. Part Number:** 8051H-900S

**2. Photograph:**



**3. Product Information:**

Description:	SteriWare® Jumbo Pipette 900mm - sterile
Use:	Single use, disposable liquid sampler
Material of Construction:	High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> <li>• Conforms to FDA 21 CFR 177.1520</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• BSE/TSE free</li> </ul> <p>A blue dye is added to one of the tips. The blue dye is certified:</p> <ul style="list-style-type: none"> <li>• FDA acceptable for food contact</li> <li>• Free of materials of animal origin</li> </ul>
Method of Construction:	Ultrasonic welding
Nominal Length	892mm
Diameter of Outer Tube:	21mm
Diameter of hole in tip:	Handle End (Blue) = 1.5mm Product End (Natural) = 3mm
Nominal Weight of Sampler:	103g
Assembly & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25kGy minimum dose)
Number of Samplers per Box:	50
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	Suitable for use in ATEX rated areas – Contact Sampling Systems for a copy of the report to confirm suitability for your application

**Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**4. Country of Origin:**

Made and packed in the UK.

**5. Repeated Use:**

This is a single-use, disposable sampler. Repeated use will lead to failure.

**6. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 10-May-19