

1. Part Number: M1044B-500

2. Photograph:



3. Product Information:

Description:

Metal Detectable Jug – 500ml - Blue

Material of Construction:

Polypropylene (PP)

- Metal detectable grade
- Only virgin material is used
- PP conforms to FDA 21 CFR 177.1520
- Blue colourant conforms to FDA 21 CFR 178.3297
- Filler conforms to FDA 21 CFR 186.1374
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

Method of Construction:

Injection moulded as a single piece

Nominal Volume:

500ml

Brimful Volume:

635ml

Top Diameter (excl. spout):

90mm

Max Width (inc. spout and handle)

135mm

Base Diameter:

75mm

Height:

130mm



Sampling Systems Ltd.

4-5 Forge Mills Park, Coleshill, Warwickshire. B46 1JH. UK.

Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com

Nominal Weight:	95g
Moulding Environment:	Class 100,000 Cleanroom
Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Number of samplers per Pack:	4
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested
Graduation Marks:	Due to the opaque blue colour the graduation marks can not be used.

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare & Metal Detectable ranges please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Apr-19