

1. Part Number: M1047B-250

2. Photograph:



3. Product Information:

Description:	Metal Detectable Jug – 500ml - Blue
Material of Construction:	<p>Polypropylene (PP)</p> <ul style="list-style-type: none"> • Metal detectable grade • Only virgin material is used • PP conforms to FDA 21 CFR 177.1520 • Blue colourant conforms to FDA 21 CFR 178.3297 • Filler conforms to FDA 21 CFR 186.1374 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Injection moulded as a single piece
Nominal Volume:	250ml
Brimful Volume:	280ml
Height:	88mm
Lid Diameter:	76mm
Nominal Weight:	97g
Moulding Environment:	Class 100,000 Cleanroom
Packing Environment:	Class 10,000 Cleanroom

Individually Bagged? Yes (heat sealed PE/PET bag)
Number of samplers per Pack: 4
Recommended Storage Conditions: Dry and ambient temperature
Shelf Life: 5 years from the date of manufacture
ATEX: This device has not been ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement.
To be kept informed of all major changes to any product in the SteriWare & Metal Detectable ranges please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 18-Apr-19