

1. Part Number: 5065B-502

2. Photo:



Angular Beaker fitted to a Telescopic Rod.
Note: Telescopic Rod is sold separately.



Angular Beaker

3. Product Information:

Description	Angular Beaker 500ml	
Materials of Construction:	Beaker	Polypropylene - Translucent
	Rotating Elbow	PBT – Blue
Nominal Volume:	500ml	
Brimful Volume:	550ml	
Maximum Beaker Diameter:	92mm	
Overall Beaker Width:	120mm	
Nominal Beaker Height:	112mm	
Nominal Weight:	131g	
Recommended Storage Conditions:	Dry and ambient temperature	

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 25-Jan-23