

**1. Part Number:** A336-1000

**2. Photograph:**



**3. Product Information:**

Description:	Open PharmaScoop® 1000ml
Material of Construction:	316L Stainless Steel
Method of Construction:	Crevice free, all interior welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	1000ml
Overall Length:	390mm
Body Length:	270mm
Body Width:	130mm
Scoop Height:	50mm
Nominal Weight:	540g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 17-May-22

Note: All dimensions are nominal. Please contact us if any dimension is critical for your application.