


1. Part Number: 5100A-250**2. Photo:****3. Product Information:**

Description

PTFE Dipper 250ml

Product Contact Parts:

PTFE (Polytetrafluoroethylene)

- Pure, virgin PTFE
- USP Class VI compliant
- FDA conforming

Nominal Capacity:

250ml

Diameter of Cup:

70mm

Handle Length:

600mm

Overall Length:

715mm

Working Temperature Range:

-200 to +280°C

Nominal Weight:

515g

Recommended Storage Conditions:

Dry and ambient temperature

Notes:

Increase the effective reach of this sampler by using Extension Rods

- PTFE Extension Rod 600mm Part No. 5100E-600
- PTFE Extension Rod 1000mm Part No. 5100E-1000

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 1-Oct-19