

**1. Part Number:** 5021P-150

**2. Photo:**



### 3. Product Information:

Description Pump Sampler Set – PE 150ml

Set Contents:

Part No.	Description	Quantity	Material
5020A-09	Hand Pump*	1	Plastic with polymer O-ring seals
5020P-10	Sampling Tube (13m length)	1	Polyethylene
5020A-07	Tube Cutter	1	Plastic with metal blade
5021H-150	Sample Bottle 150ml	12	HDPE

\* The Hand Pump is suitable for use with tubes OD from 5mm to 8mm

Recommended Storage Conditions: Dry and ambient temperature


Dimensions of Case: 390 x 300 x 100mm

Weight of Set: 1.2 kg

- Notes:
- Once the sample is full the operator should stop operating the pump. If the user continues to operate the pump then the liquid will be drawn into the body of the pump and will leak out.
  - This pump must not be used with flammable solvents

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 5-Jul-21



**Sampling Systems Ltd.**

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.  
Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com