


	<h1>Product Specification</h1>	Specification Number	A324-200
		Revision	01

Part No. A324-200  
Description Ladle  
Nominal Volume: 200ml  
Material of Construction: 304 stainless steel



Surface Finish: Better than 0.5 microns Ra  
Overall Length: 365mm  
Bowl Diameter: 100mm  
Nominal Weight: 180g

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 22-Apr-13
----------------------------	----------------------------------------------------------------------------------------------------	----------------------------------------------------------------