


1. Part Number: 7000S-600

2. Photograph:



3. Product Information:

Description	Viscous Sampler SS 600mm
Material of Construction:	316L stainless steel Silicone O-ring (FDA acceptable CFR 177.2600)
Method of Construction:	All parts screw together
Surface Finish:	Better than 1.0 microns Ra
Diameter of Outer Tube:	25mm
Diameter of hole in tip:	10mm
Max. Length:	1245mm
Min Length:	685mm
Nominal Weight:	1.4 kg
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 1-Mar-22