



	<h1 style="margin: 0;">Sampler Specification</h1>	Specification Number	7000P-1000
		Revision	01

Part No. 7000P-1000  
Description Viscous Sampler  
Material of Construction: Polypropylene (for product contact parts)



Diameter of Outer Tube: 25mm  
Max. Diameter of Tip: 33mm  
Overall Length: 1170mm  
Nominal Weight: 485g

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 29-Apr-13
----------------------------	--	--