

**1. Part Number:** 8079A-450

**2. Product Information:**

Description: PureBag Stand-Up 450ml - Sterile  
Use: Single use, disposable sampling bags  
Material of Construction: Biax Nylon/Light Density Polyethylene (LDPE) - virgin  
• Conforms to FDA 21 CFR 177.1520  
• BSE/TSE free  
Closure Type: 2 round wires with polypropylene safety tabs  
Bag Dimensions: 102 x 229mm  
Thickness of Bag: 102 microns  
Nominal Bag Volume: 400ml  
Method of Sterilisation: Gamma Irradiation  
Validation tests based on ISO 11137-1 and 11137-2 standards are performed in our facilities to confirm the efficiency of our processes.  
Packaging: 1000 bags per case (2 x 500)  
Recommended Storage Conditions: Dry and ambient temperature  
Shelf Life: 5 years  
ATEX: This product has not been ATEX tested

**3. Country of Origin:** Made and packed in Canada.

**4. Repeated Use:** This is a single-use, disposable product and should not be re-used.

**5. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 10-Nov-20

Note: The information given in this document has been given to us by our suppliers and is passed on in good faith, we believe the information to be correct.