

1. Part Number: 1360A-601

2. Photo:



3. Product Information:

Description: AccuPod® Sample Thief 600mm, 1 port
 Material of Construction: 316L stainless steel
 Method of Construction: Welds ground and polished
 Surface Finish: Better than 0.8 microns Ra (tube exteriors)
 Overall Length: 600mm
 Outer Tube Diameter: 25mm
 Number of Sampling Ports: 1
 Position of Sampling Ports:



Nominal Weight: 1.7 kg
 BSE/TSE statement: All polishing compounds used in the manufacture of this sampler are of vegetable origin
 Recommended Storage Conditions: Dry and ambient temperature
 First Use: This sampler should be thoroughly cleaned before being used for the first time
 Note: This sampler is designed for use with AccuPod® inserts. AccuPod® inserts must be ordered separately. Contact Sampling Systems for details.

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 9-Mar-20