

1. Part Number: 1360A-1501

2. Photo:



3. Product Information:

Description:	AccuPod® Sample Thief 1500mm, 1 port
Material of Construction:	316L stainless steel
Method of Construction:	Welds ground and polished
Surface Finish:	Better than 0.8 microns Ra (tube exteriors)
Overall Length:	1500mm
Outer Tube Diameter:	25mm
Number of Sampling Ports:	1
Position of Sampling Ports:	



Nominal Weight:	3.8 Kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature
First Use:	This sampler should be thoroughly cleaned before being used for the first time
Note:	This sampler is designed for use with AccuPod® inserts. AccuPod® inserts must be ordered separately. Contact Sampling Systems for details.

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 9-Mar-20