

**1. Part Number:** 1360A-2001

**2. Photo:**



**3. Product Information:**

Description: AccuPod® Sample Thief 2000mm, 1 port  
 Material of Construction: 316L stainless steel  
 Method of Construction: Welds ground and polished  
 Surface Finish: Better than 0.8 microns Ra (tube exteriors)  
 Overall Length: 2000mm  
 Outer Tube Diameter: 25mm  
 Number of Sampling Ports: 1  
 Position of Sampling Ports:



Nominal Weight: 4.9 Kg  
 BSE/TSE statement: All polishing compounds used in the manufacture of this sampler are of vegetable origin  
 Recommended Storage Conditions: Dry and ambient temperature  
 First Use: This sampler should be thoroughly cleaned before being used for the first time  
 Note: This sampler is designed for use with AccuPod® inserts. AccuPod® inserts must be ordered separately. Contact Sampling Systems for details.

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 9-Mar-20