


**1. Part Number:** 1360A-1003

**2. Photo:**



**3. Product Information:**

Description: AccuPod® Sample Thief 1000mm, 3 ports  
 Material of Construction: 316L stainless steel  
 Method of Construction: Welds ground and polished  
 Surface Finish: Better than 0.8 microns Ra (tube exteriors)  
 Overall Length: 1000mm  
 Outer Tube Diameter: 25mm  
 Number of Sampling Ports: 3  
 Position of Sampling Ports:



Nominal Weight: 2.5 kg  
 BSE/TSE statement: All polishing compounds used in the manufacture of this sampler are of vegetable origin  
 Recommended Storage Conditions: Dry and ambient temperature  
 First Use: This sampler should be thoroughly cleaned before being used for the first time  
 Note: This sampler is designed for use with AccuPod® inserts. AccuPod® inserts must be ordered separately. Contact Sampling Systems for details.

**4. Approval:**

Signed:



Name: Kuldeep Flora  
 Position: Quality Manager  
 Date: 9-Mar-20