

**1. Part Number:** 9370H-100S

**2. Photograph:**



### 3. Product Information:

|                           |                                                                                                                                                                                                                                                                                                                           |
|---------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description:              | SteriWare® ULTRA Liquid Funnel 100mm - Sterile                                                                                                                                                                                                                                                                            |
| Use:                      | Double bagged, single use, liquid funnel                                                                                                                                                                                                                                                                                  |
| Material of Construction: | 100% High Density Polyethylene (HDPE) <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.1520</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• Conforms to USP Class VI</li> <li>• Conforms to USP 661</li> <li>• Certified BSE/TSE free</li> </ul> |
| Method of Construction:   | Single piece, injection moulded                                                                                                                                                                                                                                                                                           |
| Max. Width:               | 117mm                                                                                                                                                                                                                                                                                                                     |
| ID at Top of Funnel:      | 100mm                                                                                                                                                                                                                                                                                                                     |
| Overall Height:           | 103mm                                                                                                                                                                                                                                                                                                                     |
| OD at Bottom of Spout:    | 11mm                                                                                                                                                                                                                                                                                                                      |
| ID at Bottom of Spout:    | 8mm                                                                                                                                                                                                                                                                                                                       |
| OD at Top of Spout:       | 17mm                                                                                                                                                                                                                                                                                                                      |
| Nominal Weight:           | 16g                                                                                                                                                                                                                                                                                                                       |
| Moulding Environment:     | Class 100,000 Cleanroom                                                                                                                                                                                                                                                                                                   |
| Packing Environment:      | Class 10,000 Cleanroom                                                                                                                                                                                                                                                                                                    |

|                                 |                                                                                        |
|---------------------------------|----------------------------------------------------------------------------------------|
| Bags:                           | - Double Bagged, heat-sealed PE/PET bag.<br>- Inner bag printed with batch information |
| Outer Packaging:                | Cardboard box                                                                          |
| Box Liner:                      | LDPE bag                                                                               |
| Quantity per Box:               | 20                                                                                     |
| Method of Sterilisation:        | Gamma Irradiation (25-45kGy)                                                           |
| Recommended Storage Conditions: | Dry and ambient temperature                                                            |
| Shelf Life:                     | 5 years from the date of manufacture                                                   |
| ATEX:                           | This device has not been ATEX tested                                                   |

#### 4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

#### 5. Country of Origin:

Made and packed in the UK.

#### 6. Repeated Use:

This is a single-use, disposable device. Repeated use will lead to failure.

#### 7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 6-Dec-21