

1. Part Number: 8008P-2000S

2. Photograph:



3. Product Information:

Description:	SteriWare® Graduated Beaker 2000ml - Sterile
Use:	Single use, disposable beaker
Material of Construction:	Polypropylene (PP) - virgin <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Nominal Volume:	2000ml
Brimful Volume:	2400ml
Max. Diameter (inc. spout):	165mm
Top diameter (excl. spout)	146mm
Bottom Diameter:	121mm
Overall Height:	195mm
Graduation Marks:	From 100ml to 2000ml (every 100ml)
Nominal Weight:	161g
Moulding Environment:	Class 100,000 Cleanroom

Primary Packing Environment: Class 10,000 Cleanroom
Individually Bagged? Yes (heat sealed PE/PET bag)
Method of Sterilisation: Gamma Irradiation (25 to 45 kGy)
Quantity per Box: 15
Recommended Storage Conditions: Dry and ambient temperature
Shelf Life: 5 years from the date of manufacture
ATEX: This device has not been ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement. To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable device. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 6-Mar-23