

# Protease Inhibitor Mix FY

[Pack Sizes](#)

[Certificates of Analysis \(Lot.-no. - release date\)](#)

[PDF Documents](#)

[Related Products](#)

Protease Inhibitor Mix FY is an effective protection against proteolytic degradation in the preparation of fungal and yeast extracts. The mixture of four protease inhibitors has a broad range of activity against aspartate-, cysteine-, serine-, and metallo proteases. Protease Inhibitor Mix FY contains AEBSF, E-64, Pepstatin A, and 1,10-Phenanthroline. Supplied as a kit: each vial is provided with an extra vial of 1 ml DMSO. The content of 1 vial dissolved in 1 ml DMSO gives a 100-fold concentrate, sufficient for 99 ml extract.

DANGER

Hazard Statements

[H301 - H410](#)

Precautions

[P273](#)

Reaction

[P301 +P310 - P330 - P391](#)

Storage

[P405](#)

GGVSE/ADR: 6.1 III UN2811 • IATA: 6.1 III UN2811 • EINECS: 200-664-3 • WGK: 1 •

HS: 38229000

Storage Temperature: -15 °C to -25 °C

