

Référence produit: 127263-2DG



— Informations principales

Description du produit:	Cryo.s, pp, steril, 4,0 ml, datamatrix
Stérile:	sterile
Conditionnement:	300
Sous conditionnement :	50

— Emballage

Poids emballage:	1,37 kg
Dimension emballage:	323 x 243 x 127 mm
Conditionnement:	300
Sous conditionnement :	50
PAL:	50400

— Informations détaillées

Volume de remplissage (ml)\n:	0.00
-------------------------------	------

Cryo.s Freezing Tubes, 4 ml – with Datamatrix Coding

- Free of detectable DNase, RNase, human DNA
- Non-Pyrogenic, non-cytotoxic
- Endotoxin-free
- CE-marked
- Made of high-grade polypropylene which guarantees high thermal resistance in a temperature range from -196°C to +121°C
- As a safety precaution Cryo.s™ should be stored in a freezer or exclusively in the gas phase over the liquid nitrogen
- Suitable for shipment by air
- Laser-written ECC200 Datamatrix
- Strong resistance to mechanical and chemical damage
- High data capacity
- Only Datamatrices with readability grades A and B pass quality control
- Round bottom
- Working volume: up to 4.0 ml
- Diameter: 12.4 mm; total height incl. lid: 83 mm
- With starfoot base and external thread
- With screw cap in natural
- Sterile
- Description 2D Barcode: 14x14 Datamatrix, sequence: G700000001, G700000002 ...

Quantity per bag/case: 50/300

DOMINIQUE DUTSCHER SAS