


EUKITT® *classic*

The gold-standard of mounting media for cover slipping machines



RELIABLE

GOLD-
STANDARD

EUKITT® classic

Mounting Media

- Fast drying for rapid results
- Neutral and colourless
- Well proven formulation
- Long-term stability

EUKITT® classic has proven long term storage qualities. Slides exhibit minimal discolouration or drying back. Preparations mounted with EUKITT® classic lead to superior optical quality for reliable long-term filing. Ideal viscosity formulation for automatic cover slipping instruments. No bubbles or strikes due to ideal vapor pressure.

Supplied in a wide neck, unbreakable aluminium bottle. The product has a shelf life of 3 years.



Technical specifications:

Refractive index $n[20\text{ °C}]$: 1.49

Viscosity: 500 – 750 mPa*s

Drying time: 20 minutes

No yellowing after long storage

No brittle fracture

No self-fluorescence under UV light

EUKITT® classic is available in the following packaging units:

100 ml Bottle	250 ml Bottle	500 ml Bottle
6 x 100 ml	10 x 250 ml	6 x 500 ml
10 x 100 ml	25 x 250 ml	10 x 500 ml
25 x 100 ml		25 x 500 ml

For further information please visit us:

www.eukitt.net



is now a division of
ORSA-tec®

Contact

ORSAtec GmbH
Max-Fischer-Straße 11
86399 Bobingen

Tel.: +49 (0) 761 81077
Fax +49 (0) 761 892535
info@eukitt.de