


PDS No. 1242xx	PRODUCT DATA SHEET				Page 1 of 1
Revision 09	Cryo.s™, 5 ml, Round Bottom			CE	
	Item-No. 1242xx				
Valid for Item-No.:	124261 (sterile)	124263 (sterile)	124273 (sterile)	124274 (sterile)	124275 (sterile)
	124276 (sterile)				

1.	Description / Specification	
1.1	Description	Cryo.s™ with different coloured screw caps (caps have a silicone gasket), internal thread, round bottom and additional cap inserts for labelling 124261: without writing area 124263, -273, -274, -275, -276: with writing area and graduation Additional inserts, Cat.-No. 304 171 (50 pieces per bag)
1.2	Dimensions	See Customer Drawing
1.3	Volume	Working volume: 4,5 – 5,0 ml
1.4	Material / Resin	<u>Tube</u> : PP (Polypropylene), free of heavy metal <u>Screw cap</u> : PP (Polypropylene), free of heavy metal <u>Gasket</u> : Silicone <u>Small cap insert</u> : PS (Polystyrene), free of heavy metal
1.5	Colour	<u>Tube</u> : natural, writing area and graduation: white <u>Screw cap</u> : 124261, -263: natural 124273: red 124274: blue 124275: green 124276: yellow <u>Gasket</u> : translucent <u>Small cap insert</u> : white
1.6	Sterilization	SAL 10 ⁻⁶
1.7	Quality Control	<u>Raw Material-Control</u> : physical testing <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2.	Temperature range	-196°C to +100°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	5800 x g: swinging-bucket rotor 34000 x g: fixed-angle rotor
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	5 years
2.7	Other Information	Recommended for storage in the freezer or above liquid nitrogen

3.	Packaging	
3.1	Pieces / Bag	50
3.2	Pieces / Box	300
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	<u>Additional in each box</u> : Safety precaution 2 bags with 50 cap inserts per bag (100 pieces) Certificate of Quality to download

4.	Other Information
4.1	Laboratory use only

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 08	Date 16 December 2019	Date 17 December 2019	Date 17 December 2019	
Date 19.06.2019	Name S. Kaelberer	Name Dr. S. Mühlfriedel	Name A. Illig	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.