

Référence produit: 126263-2DG



#### — Informations principales

Description du produit:	CRYO.S PP, STERIL, 2 ML 12,4/47 MM,DA
Stérile:	sterile
Conditionnement:	500
Sous conditionnement :	100

#### — Emballage

Poids emballage:	1,24 kg
Dimension emballage:	323 x 243 x 127 mm
Conditionnement:	500
Sous conditionnement :	100
PAL:	84000

#### — Informations détaillées

Volume de remplissage (ml)\n:	0.00
-------------------------------	------

## Cryo.s Freezing Tubes, 2 ml – with Datamatrix Coding

- Free of detectable DNase, RNase, human DNA
  - Non-Pyrogenic, non-cytotoxic
  - Endotoxin-free
  - CE-marked
  - Made of high-grade polypropylene which guarantees high thermal resistance in a temperature range from -196°C to +121°C
  - As a safety precaution Cryo.s™ should be stored in a freezer or exclusively in the gas phase over the liquid nitrogen
  - Suitable for shipment by air
  - Laser-written ECC200 Datamatrix
  - Strong resistance to mechanical and chemical damage
  - High data capacity
  - Only Datamatrices with readability grades A and B pass quality control
  - Round bottom
  - Working volume: up to 2.2 ml
  - Diameter: 12.4 mm; total height incl. lid: 47 mm
  - With starfoot base and external thread
  - With screw cap in natural
  - Sterile
  - Description 2D Barcode: 14x14 Datamatrix, sequence: G600000001, G600000002 ...
- 

DOMINIQUE DUTSCHER SAS