



**PRODUCT CODE: 126352**

**Amido Black 10B solution for the determination of protein content in milk for analysis**

according to NF V 04-216

---

TARIC 3822 00 00 00

**PHYSICAL DATA:** liquid, Deep blue, Miscible with water • D 20/4 1,007 •

**BIBLIOGRAPHY:**

**HAZARDOUS:**

**WEIGHT/VOLUME INFORMATION:** 1l~1,007 kg 1kg~0,993 l

**OBSERVATIONS:** Storage between 15 and 25°C • Storage away from direct light. •

**SPECIFICATIONS:**

pH at 20°C 2,30-2,50

ABS  
at  $\lambda$  620 nm 1% sol. in H<sub>2</sub>O 0,695-0,735