


CATALOG NO.....:	4300-0250	
SIZE (NOM. VOL. ml) ...:	250	
NOM. BRIM CAP (ml) .:	364	
NOM. WEIGHT (gms) ...:	71	
A (NECK I.D.).....:	INCH	MM
B (HT. W/ CLOSURE).....:	0.95	24.1
C (HT. W/O CLOSURE).....:	12.86	326.6
D (O.D. AT BULB AREA) .:	12.71	322.8
E (TUBING I.D.).....:	3.02	76.7
W.T. (MIN WALL THK) .:	0.28	7.1
	0.035	0.89

MATERIALS	
CLOSURE	PP
FUNNEL BODY	PP
STOPCOCK ASSEMBLY	TEFLON TFE
TIP	TEFLON FEP

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

**ThermoFisher**  
S C I E N T I F I C

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

NAME	TL	GROUP	GRP
DATE	GP		5/11/09
DESCRIPTION:	PSS 00-4300-0250		

REV. 0409

SEPARATORY FUNNEL

**CONFIDENTIAL**

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.